

# DISINFECTION OF AMBULANCES WITH NOCOSPRAY®/NOCOLYSE®

## POSSIBILITY 1:


Use the Nocospray® as in a standard room : place the device directly in the vehicle



### POSSIBILITY 2:

Use a Nocospray® with a wall-crossing system

⇒ Enables to treat several cars with only one Nocospray quickly



### POSSIBILITY 3:

Use a wall-mounted Nocospray®

⇒ fixed on the wall with the bottle in the same room as the diffusion



### POSSIBILITY 4:

Use a wall-mounted Nocospray® with wall-crossing system

- ⇒ fixed on the wall with the bottle in the different room than the diffusion (to disinfect the vehicle without entering it)

